AUTOCLAVE LOG SHEET

LOCATION (Building/Room)						AUTOCLAVE (Make/Model)			
DATE	LOAD CONTENTS*		TIME (Minutes)	PRESSURE (psi)	MAX TEMP REACHED	TAPE RESULT (Pass/Fail)	BIOLOGICAL INDICATOR USED? (Y/N)	OPERATOR	COMMENTS

- Loads containing biohazardous waste must be autoclaved at ≥ 121°C with a sterilization time of at least 30 minutes.
- Biological Indicators must be performed at least once monthly.

